

## Testing different batches of collodion

Medium: 13x18 cm Black Acrylics

Sensitizing Time: 3:30min

Developer: Ferrosulfate + Glaecial Acid

Fixer: Plain Hypo

Plates are unvarnished and not post processed except for setting the black point to the naked black of the plates.

1) KI + KBr Collodion from 25/11/2012, 10 seconds exposure, 17 seconds development

2) Old Workhorse (MAMUT Photo, Czech) from 26/11/2012, 10 seconds exposure, 15 seconds development

3) New Guy Negative (MAMUT Photo, Czech) from 26/11/2012, 10 seconds exposure, 17 seconds development

4) KI + KBr Collodion from August 2012, 15 seconds exposure, 18 seconds development

5) Lea Landscape #7 (MAMUT Photo, Czech) from 26/11/2012, 10 seconds exposure, 15 seconds development

I/O  
AMBROTYPE &  
PHOTOGRAPHER

